trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

ENGLISH	3
DEUTSCH	18
FRANÇAIS	32
NEDERLANDS	47
ITALIANO	62
NORSK	76
SUOMI.....	91
DANSK.....	106
SVENSKA	120
ESPAÑOL	135
TÜRKÇE	149

trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

Read this addendum together with the trophon®2 user manual before operating device to determine the correct procedures.

The Table of Contents in this addendum lists updates to the trophon®2 User Manual.

For further information, contact your customer service representative or visit the Nanosonics website.

© 2026 Nanosonics Limited. All rights reserved.

www.nanosonics.com


Table of Contents

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS.....	6
A1.1 Labels and Symbols.....	6
A2.1 Indications for Use	9
Part B – SETUP	10
B2.3 Initial Setup	10
B2.4 Warm-up Cycle	10
B2.6 Basic Settings Network	10
Part E – MAINTENANCE AND ROUTINE CARE	12
SECTION E4: Cleaning and Disinfecting.....	12
Part H – CYBERSECURITY	13
APPENDIX 1: trophon2 device Technical Specifications	16
GLOSSARY	17

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS

The following section of the labels and symbols table has been updated to add a new feature relating to DICOM in Touch Screen Symbols.

A1.1 Labels and Symbols

	Caution		Warning
	Consult Instructions For Use		Corrosive
	Environmental Conditions: trophon Storage and Transport Conditions: Temperature range: -20 °C to +60 °C / -4 °F to +140 °F		Single Use Only
	Fragile / Handle With Care		UN 2014 – Hydrogen Peroxide
	Do not disassemble		Dangerous Voltage
	Separate collection for electrical and electronic equipment.		Keep Dry
	Keep Out of Direct Sunlight		This Way Up
	Batch Number		Product Number
	Serial Number		Expires (year and month)
	Legal Manufacturer		Date of Manufacture
	Oxidizer – 5.1		Corrosive – 8
	Warning: Hot Surface		Warning: Moving parts, do not touch mechanism
	Cannot be transported by air freight		Wear Gloves

	Environmental Conditions: Operating temperature range for the trophon device: 17 °C to 27 °C / 62.6 °F to 80.6 °F		AcuTrace® RFID Zone
	EU Importer		European Authorised Representative
	Conforming to Medical Devices Regulation (EU)2017/745 (MDR)		Conforming to RoHS 3 (EU 2015/863)
	Medical Device		Conforming to UK MDR (SI 2002 No 618)
Touch Screen Symbols			
	Start up from sleep		Cycle Start
	Connected to network		Disconnected from network
	Connected to the cloud		Disconnected from the cloud
	trophon software download		Menu
	No Data Transmission to PACS Server		
Integrated Probe Positioning Guide			

The following section of the indications for use has been updated relating to the Minimum Operational Cycle Time.

A2.1 Indications for Use

The trophon2 device is intended for the high-level disinfection (HLD) of non-lumened, reusable, transiently invasive and non-invasive medical instruments/devices* e.g. devices that are intended for use for imaging, diagnostic, ablation, coagulation and their accessories.

The trophon2 system consists of a multiple-use device, combined with a single-use disinfectant “trophon NanoNebulant”, delivered from a multi-dose cartridge.

The trophon2 device is suitable for use in general hospital and health care facilities by trained personnel.

The trophon NanoNebulant should be used with the following contact conditions:

Minimum Operational Cycle Time: 4 minutes

Minimum Concentration: 31.5%

Minimum Disinfectant Dose: 1.0 g

Minimum Chamber Temperature: 56 °C

The trophon2 device is NOT intended to reprocess single-use probes or instruments, or to pre-clean medical instruments.

Chemical indicator use is required with every HLD cycle. Only the trophon Chemical Indicator product is the approved chemical indicator for use with trophon2 device.

* The terms “ultrasound probe” and “probe” in the User Manual refer to approved medical instruments/devices.

Part B – SETUP

The following Installation Guide section of the user manual has been added for a new feature related to DICOM.

B2.3 Initial Setup

To access the trophon2 device software version and build number:

Select **Menu -> Information -> Device**

B2.4 Warm-up Cycle

1. The warm-up cycle prepares the trophon2 device for operation and will begin automatically when the device is powered on.
2. The screen message will indicate when the trophon2 device is ready for use. Follow the onscreen instructions.
3. It is possible to load a probe ready for disinfection during the warm-up cycle, between cycles. Only load the probe when the 'Load Probe' feature is displayed on the screen. When the available, select 'Load probe' to open the chamber door and follow the onscreen instructions. The cycle shall begin once the chamber reaches the appropriate temperature.

NOTE: The HLD cycle cannot begin until the required minimum chamber temperature is met. (see section A2.1 Indications for Use).

B2.6 Basic Settings Network

For DICOM user, follow additional instructions as below:

- **DICOM**

To view DICOM network settings, select:

Menu → Settings -> Network -> DICOM

When a trophon2 Plus Software Upgrade device is not connected to a DICOM PACS server, it is indicated at the top right of the idle screens, by the following symbol:

To configure DICOM server list, select:

Menu → Settings -> Network -> DICOM -> Servers

To configure DICOM server settings, select:

Menu → Settings -> Network -> DICOM -> Services

To check subscription and other DICOM settings, select:

Menu → Settings -> Network -> DICOM -> Setting

To clear DICOM cache records and restore default DICOM configuration, select:

Menu → Settings -> Network -> DICOM -> Maintenance

To enable or disable DICOM traceability to the Nanosonics cloud service, select:

Menu → Settings -> Network -> DICOM -> Traceability Data

NOTE: Mandatory entry fields are indicated where an * is detailed on the screen.

Part E – MAINTENANCE AND ROUTINE CARE

SECTION E4: Cleaning and Disinfecting

- Do NOT submerge the trophon2 device, or pour liquids over it.
- Keep the trophon2 device level and upright at all times.
- Keep the power socket completely dry (see Figure 4).

Ensure appropriate PPE, including clean gloves, is worn when cleaning and disinfecting the trophon2 device.

For cleaning, Nanosonics recommends that the surfaces of the trophon2 device be thoroughly cleaned in accordance with organisational policies and following the detection of visible contamination. When cool, thoroughly wipe the chamber including the silicone seals, IPP and cable clamp and outside surfaces of the trophon2 device including the screen, handle and underneath the cable tray (refer to section B1.2 Cable Tray for tray removal) with a cloth or wipe moistened with a mild general purpose detergent, or soap and water solution, such as dish soap or general cleaners, until all surfaces are visibly clean. The cleaning process should be repeated if the trophon2 device is not visually clean. Refer to the product's IFU for general instructions and directions for use.

For disinfecting, Nanosonics recommends that the surfaces of the trophon2 device be disinfected in accordance with organisational policies. Wipe all accessible surfaces of the trophon2 device with an Isopropanol or Quaternary Ammonium (Quat) wipe. Refer to the wipe product's IFU for general instructions and directions for use. Ensure surfaces remain wet for the full contact time for disinfection according to the wipe manufacturer's IFU and use additional wipes if required.

Part H – CYBERSECURITY

The trophon2 device has built-in cybersecurity features designed to help protect both the device and your hospital network. To maintain the highest level of protection, ensure the device is updated to the latest software version.


It is recommended to connect the trophon2 Device to the hospital network via the ethernet, WiFi or Cellular to ensure timely software updates that provide the latest features and security enhancements.

Consult your IT department to network your device at time of installation. If required, contact your customer service representative to request a copy of the Manufacture Disclosure Statement for Medical Device Security (MDS2).

trophon2 Cybersecurity Architecture

The trophon2 Cyber Architecture (Figure 1) gives a pictorial view of the connections made by and to the Trophon2 device in a network environment.

Figure 1: Cybersecurity Architecture

trophon2 Software Updates

trophon2 device software update methods:

- Software updates can be applied either by a certified service technician on-site, via self-serve tool, or remotely when the device is connected to the Nanosonics cloud service.
- When a new update is available, the device notifies the user via the touchscreen. Users can choose whether to proceed with the update.

System Monitoring and Log Management

The trophon2 device automatically records operational activity and any abnormal events in internal logs. These logs are accessible only to Nanosonics representatives using specialized tools.

Consult Nanosonics Service if you notice any inconsistencies or strange behaviour from the trophon2 Device.

trophon2 Device protective features

Communications between the device and approved service applications are encrypted and restricted to essential functions. Only connections with valid security certificates are accepted. All other network traffic is blocked.

Backup and Restore features / restoring authenticated user

Your trophon2 Device configuration is verified along with the device software to create an intrinsic set of parameters used to operate the device. Since the device configuration is part of the software, there is no need to backup, restore or change this configuration.

Log Files

Any abnormal events will be recorded in the device logs. These are not user accessible but may be accessed by a Nanosonics representative using appropriate software on an external device.

trophon2 Device end of life

The trophon2 Device software is prepared by Nanosonics and does not depend upon external components for operation. As such, Nanosonics can ensure continued support for upgrades of the device software into the future.

Secure decommission of the trophon2 Device

There is no patient information stored on the trophon2 Device. Decommissioning of the device is acceptable without erasure of data.

Additional Information

A copy of Software Bill of Materials (SBOM), Manufacture Disclosure Statement for Medical Device Security (MDS2), and other software license details can be provided on request by contacting Nanosonics Customer Service at customerservice@nanosonics.us.

trophon2 Device Software Licenses

The trophon2 device software uses statically linked components that are under MIT, GPL-2.0, LGPL 2.1, BSD 3-clause, EPL 1.0, GPL-3.0 Apache-2.0 licenses.

- MIT license [1]: Yocto Project-Dunfell, libudev

- GPL 2.0 license [2]: Linux Kernel, U-Boot

- Apache-2.0 license [3]: Libraries used in Java and Qt backend, eg; spring, grpc, netty, httpclient, dbus, quartz, log4j, derby etc

- BSD 3-clause [4]: protobuf-java, libupb
- LGPL 2.1 [5]: hibernate, libz
- EPL 1.0 [6]: junit
- GPL 3.0 [7]: libgpr

Please see below for links to the referenced license notes:

- [1] <https://opensource.org/license/mit>
- [2] <https://www.gnu.org/licenses/old-licenses/gpl-2.0.txt>
- [3] <https://www.apache.org/licenses/LICENSE-2.0.txt>
- [4] <https://opensource.org/license/bsd-3-clause>
- [5] <https://www.gnu.org/licenses/old-licenses/lgpl-2.1.en.html>
- [6] <https://www.eclipse.org/legal/epl/epl-v10.html>
- [7] <https://www.gnu.org/licenses/gpl-3.0.en.html>

APPENDIX 1: trophon2 device Technical Specifications

Corrected typographical error in electromagnetic compliance section for product class.

Electrical Specification	Rated mains input voltage: 230 V AC Rated mains input current: 6 A, 50/60 Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using the power cable supplied with the trophon3 device.
Data port	Ethernet connector RJ45 USB Port: Type A
Environmental Specification	Operating temperature range: 17 °C to 27 °C / 62.6 °F to 80.6 °F
Storage and Transport Conditions	Temperature range: -20 °C to +60 °C / -4 °F to +140 °F
Physical Characteristics	Weight of trophon2 device: Unpacked: 22 kg (48.5 lbs) Packed: 25 kg (55 lbs) Dimensions of trophon2 device: 535 mm high × 360 mm wide × 317 mm depth (21 in high × 14.2 in wide × 12.5 in depth)
Electromagnetic Compliance	The trophon2 device has been tested and found to comply with the limits for emission (electromagnetic Interference) pursuant to EN61326-1:2013 (CISPR 11 Group 1 Class A limits)

GLOSSARY

PACS

Picture Archiving and Communication System

DICOM

Digital Imaging and Communications in Medicine

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
**Manufactured by:**

Nanosonics Limited
3-5 Thomas Holt Drive,
Macquarie Park NSW 2113,
Australia
+61 2 8063 1600
customerservice@nanosonics.com.au
www.nanosonics.com.au

**Authorised European Representative/
Importer**

Nanosonics Europe GmbH
Poppenbütteler Bogen 66
22399 Hamburg
Germany
+49 40 46856885
customerservice@nanosonics.eu
www.nanosonics.eu

UK Responsible Person:

Nanosonics UK Limited
Ground Floor at The Forum
Unit C1 & C2, Hercules Business Park,
Bird Hall Lane, Stockport, SK3 0UX, UK
+44 (0) 161 686 3030
customerservice@nanosonics.co.uk
www.nanosonics.co.uk

Conforming to WEEE Directive 2012/19/EU under Article 7 Recovery

www.nanosonics.com