trophon[®]2

User Manual

trophon[®]2

User Manual

Read this manual before operating the trophon®2 to determine the correct procedures.

For further information, contact your customer service representative or visit the Nanosonics website.

All technical specifications and system approvals are listed in Appendix 1.

©2018 Nanosonics Limited. All rights reserved.

The content of this manual is correct at the time of product purchase.

trophon®2, NanoNebulant® and Sonex-HL® are registered trademarks of Nanosonics Limited.

AcuTrace™ is a registered trademark of Nanosonics Limited in Australia and other countries.

NanoNebulant is the product name of the trophon2 disinfectant used in all regions where trophon2 is available for sale, with the exception of US and Canada.

Sonex-HL is the product name of the trophon2 disinfectant in the US and Canada.

Nanosonics' technologies are protected by a combination of patents, trademarks, and exclusive proprietary rights: <http://www.nanosonics.com.au/nanosonics-patent>

Your trophon2 representative is:

Attach Business Card or information sticker/stamp here.

 ⁰¹⁹⁷ Conforming to MDD 93/42/EEC; certified by Notified Body TUV Rheinland

Conforming to RoHS Directive 2011/65/EC

 Conforming to WEEE Directive 2012/19/EU under Article 7 Recovery

www.nanosonics.eu • www.nanosonics.com.au

Table of Contents


Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS	5
SECTION A1: Important Labels, Symbols and Warnings	5
A1.1 Labels and Symbols.....	5
A1.2: Warnings.....	6
SECTION A2: Introduction to the trophon2.....	7
A2.1 Indications for Use	7
A2.2 Disinfection Process	7
A2.3 Validated Probes, Disinfectants and Chemical Indicators.....	7
A2.4 Training.....	8
A2.5 Environment and User Profile.....	8
SECTION A3: Instructions.....	8
Part B – SET UP.....	9
SECTION B1: trophon2 Overview.....	9
B1.1 trophon2 Features	9
B1.2 Cable Tray	10
SECTION B2: Installation Guide	11
B2.1 Positioning your trophon2.....	11
B2.2 Powering On.....	12
B2.3 Initial Setup.....	12
B2.4 Warm up Cycle	12
B2.5 trophon2 Touch Screen	12
B2.6 Basic Settings	12
B2.7 AcuTrace™	12
B2.8 AcuTrace™ Settings	14
SECTION B3: trophon AcuTrace PLUS.....	14
B3.1 Activation	14
B3.2 Network Parameters Setup	14
PART C – OPERATION	14
SECTION C1: Loading the Disinfectant Cartridge	14
SECTION C2: Logging the trophon Chemical Indicators.....	14
SECTION C3: Routine High Level Disinfection Cycle	15
C3.1 Preparing the Probe	15
C3.2 Inserting the Chemical Indicator.....	15
C3.3 Positioning the Probe	15
C3.4 Closing the Chamber Door	17
C3.5 Disinfecting the Probe	17
C3.6 Removing the Probe.....	18
C3.7 Sleep Mode	18
PART D – RECORDS.....	18
SECTION D1: Record Options	18
PART E –MAINTENANCE AND ROUTINE CARE	19
SECTION E1: Preventative Maintenance Service	19
SECTION E2: Purge Cycle.....	19
E2.1 When to Run a Purge Cycle	19
E2.2 How to Initiate a Purge Cycle	20
SECTION E3: Regular Cleaning.....	20
SECTION E4: Transporting the trophon2	20
SECTION E5: Disposal of trophon2	20
PART F –TROUBLESHOOTING.....	20
SECTION F1: Incomplete or Failed Cycles	20


F1.1 Mains Power Failure	20
F1.2 Cycle Fault	21
SECTION F2: Manual Door Lock Override.....	21
SECTION F3: Diagnostics	22
 PART G – SERVICE REQUIREMENT AND WARRANTY PROVISION	24
 APPENDIX 1: TROPHON2 TECHNICAL SPECIFICATIONS N05000-2	24
 APPENDIX 2: PRODUCT WARRANTY TERMS AND CONDITIONS ..	25
 APPENDIX 3: MicroDoc LICENCE AGREEMENT	26
 APPENDIX 4: TROPHON2 ACCESSORIES*	27
 GLOSSARY	28


Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS

SECTION A1: Important Labels, Symbols and Warnings

A1.1 Labels and Symbols

	Caution		Warning
	Consult instructions for use		Corrosive
	Environmental Conditions: trophon2 Storage and Transport Conditions: Temperature range: -20°C to +60°C / -4°F to +140°F		Single Use Only
	Fragile / Handle With Care		UN 2014 – Hydrogen Peroxide
	Do not disassemble		Dangerous Voltage
	Separate collection for electrical and electronic equipment		Keep Dry
	Keep Out of Direct Sunlight		This Way Up
	Batch Number		Product Number
	Serial Number		Expires (year and month)
	Legal Manufacturer		Date of Manufacture
	Oxidizer – 5.1		Corrosive – 8
	Warning: Hot Surface		Warning: Moving parts, do not touch mechanism
	Cannot be transported by air freight		Wear Gloves

	Environmental Conditions: trophon2 Operating temperature range: 17°C to 27°C / 62.6°F to 80.6°F		AcuTrace™ RFID Zone
	Start Up from sleep		Cycle Start
	Menu		

Integrated Probe Positioning Guide

A1.2: Warnings

Hot Temperatures

- Do NOT touch surfaces in the internal chamber. They can be hot and cause burns.
- To avoid probe damage, ensure the probe is correctly positioned in the chamber. (See section C3.3 for correct positioning of probe.)

Malfunctions

- Do NOT attempt to open the chamber door during a cycle, power failure or system malfunction.
- All repairs must be carried out by trained personnel.

Transporting the trophon2

Follow your facilities manual handling procedures for guidance on lifting heavy objects.

- The trophon2 weighs approximately:
Unpackaged 22kg (48.5lb).
Packaged 25kg (55lb).
- If your trophon2 has been used, purge the trophon2 before transportation to remove the disinfectant. (See section E2).

Electrical Device

- Use the power cable supplied with the trophon2, connect to an earthed power outlet with the correct voltage and frequency as specified on the product and in Appendix 1. Incorrect voltage can damage the product.
- The trophon2 must not be connected to the same circuit as a critical patient or life support device.
- Spilled fluid can result in electrical shock. Avoid spilling fluids on or around the trophon2. Do not immerse any parts of the trophon2 in liquid.
- Do not attempt to access the internal mechanics. This may result in electric shock.

Protective Wear and Spills

- Wear clean disposable gloves throughout the complete High Level Disinfection (HLD) process including but not limited to running the trophon2 and handling:
 - Disinfectant cartridges. Temporary bleaching and/or irritation of the skin may occur if gloves are not worn.
 - Probes before and after a HLD cycle.
 - Chemical Indicators before and after a HLD cycle.
 - Waste drawer when emptying or obtaining the Manual Door Lock key.
- Wear appropriate personal protective equipment (PPE) when managing spills.
- Never return spills to original cartridges for re-use.

SECTION A2: Introduction to the trophon2

A2.1 Indications for Use

The trophon2 is intended for the high level disinfection (HLD) of non-lumened, reusable, transiently invasive and non-invasive medical instruments/devices e.g. devices that are intended for use for imaging, diagnostic, ablation, coagulation and their accessories.

The trophon2 is designed to provide High-Level Disinfection (HLD) of validated ultrasound probes. High-Level Disinfection is achieved by surface exposure to a controlled dose of hydrogen peroxide mist delivered to a disinfection chamber containing the ultrasound probe.

The trophon2 system consists of a multiple use instrument combined with a single use disinfectant “trophon NanoNebulant”, delivered from a multi-dose cartridge.

The trophon2 is suitable for use in general hospital and health care facilities by trained personnel.

The trophon NanoNebulant should be used with the following contact conditions:

Minimum Operational Cycle Time:	7 minutes
Minimum Concentration:	31.5%
Minimum Disinfectant Dose:	1.0 g
Minimum Chamber Temperature:	56°C

The trophon2 is NOT intended to reprocess single-use devices or pre-clean medical instruments.

Chemical indicator use is required with every HLD cycle. Only the trophon Chemical Indicator product is the approved chemical indicator for use with trophon2.

** The terms “ultrasound probe” and/or “probe” in the User Manual refer to approved medical instruments.*

A2.2 Disinfection Process

At the start of the HLD cycle, the trophon2 creates an aerosol of concentrated hydrogen peroxide. This is distributed over the exposed surface of the probe providing thorough HLD of the probe’s shaft and handle. Once used, the hydrogen peroxide is converted into its constituents; oxygen and water. During a purge cycle, the converted oxygen is vented into the atmosphere and the water is collected inside the waste drawer located at the side of the trophon2 for emptying.

The contact conditions are fixed cycle parameters that the end user cannot modify.

A2.3 Validated Probes, Disinfectants and Chemical Indicators.

For details of probes that can be used in the trophon2, refer to the Validated Probe List on the Nanosonics website.

Use only trophon NanoNebulant disinfectant cartridges and trophon Chemical Indicators when high level disinfecting with the trophon2. No other disinfectant or chemical indicator has been approved for use with the trophon2.

A2.4 Training

Before setting up or using your trophon2, ensure that all users are trained in safety procedures and potential hazards, as outlined in this manual.

Visit www.nanosonicsacademy.com to complete the online training module.

A2.5 Environment and User Profile

The trophon2 is designed for use in healthcare facilities to high level disinfect ultrasound probes under the control of trained healthcare professionals.

The disinfection cartridge, Chemical Indicator, and trophon2 system are designed to be used with minimal personal protective equipment (gloves only) and in a standard workplace or clinical setting, including at the patient point-of-care. Special ventilation and other safety precautions are not required when used as per these instructions.

SECTION A3: Instructions

Read these instructions before using the trophon2:

- Safety Data Sheet (SDS) enclosed with the trophon2.
- Occupational or Workplace Health and Safety Guidelines (OH&SG, OSHA, WHS) for your institution for lifting, spills etc.
- trophon Chemical Indicator Instructions for Use (IFU).
- trophon disinfectant cartridge IFU.
- trophon accessory IFUs for any additional accessories purchased with the trophon2 (see Appendix 4).
- Probe manufacturer's instructions.

Failure to follow instructions may result in:

- Burns, bleaching, electric shock or other injury.
- High level disinfection not achieved.
- Residual disinfectant remaining on the probe, which may cause injury when removing.
- Equipment damage.

Part B – SET UP

SECTION B1: trophon2 Overview

B1.1 trophon2 Features

Front

Figure 1

Back

Figure 2

Right Side

Figure 3

Left Side


Figure 4

Chamber

Figure 5

1. AcuTrace reader.
2. User Interface (UI).
3. Chamber door handle.
4. Manual door lock opening mechanism cover.
5. Cartridge door. **Warning: Do NOT force open** (cartridge door will automatically open when required).
6. Waste drawer.
7. Power switch.
8. Power socket.
9. Ethernet port.
10. 3 x USB Ports.* **
11. Chamber door (opened).
12. Chemical Indicator holder.
13. Probe in correct position.
14. Door lock mechanism. **Warning: Do NOT put fingers into mechanism.**
15. Cable clamp.
16. Cable seal.
17. Integrated Probe Positioner (IPP).

* trophon2 does not support Wi-Fi dongles or 3G/4G dongles.

** The 3 USB ports may be used in any order.

B1.2 Cable Tray

The cable tray holds the ultrasound cable away from the chamber door and stores the cable during HLD. Remove the tray for cleaning (figure 6) and wipe with a cloth moistened with a soapy solution.

Figure 6

SECTION B2: Installation Guide



B2.1 Positioning your trophon2

 **trophon2 weighs approximately 22kg (48.5 lbs). Follow your facilities manual handling procedures for guidance on lifting heavy objects**

1. Ensure the surface is level, can support the weight and allows adequate airflow (see Appendix 1).
2. There are no specific illumination requirements to use the trophon2. Follow your facilities standard for work environment illumination guidance.
3. Ensure the area around your trophon2 is free from other equipment and clutter. Position as shown in Figure 7 to ensure access to all features.

The trophon2 can also be mounted to the wall using the trophon Wall Mount 2 or for a mobile solution, the trophon Cart. Please refer to the product IFUs for details.

 **The trophon2 should be placed at a height from the floor level to accommodate a range of user heights. Refer to Figure 8 for a guide to an ergonomically safe work zone.**

B2.2 Powering On

The trophon2 must not be connected to the same circuit as a critical patient or life support device.

1. Attach the power cable supplied to the trophon2 power socket.
2. Connect to mains power.
3. Turn on the power switch, located on the left side of the trophon2.

NOTE: To maximise the life of the disinfectant cartridge, keep your trophon2 switched on at all times, except when moving the device.

B2.3 Initial Setup

The Initial Setup launches automatically, prompting you to configure the optional settings when turning on the trophon2 for the first time. Follow the on screen instructions.

B2.4 Warm up Cycle

1. The warm up cycle prepares the trophon2 for operation and will begin automatically when the machine is powered on.
2. The screen message will indicate when the trophon2 is ready for use. Follow the on screen instructions.

B2.5 trophon2 Touch Screen

The trophon2 is operated using the touch screen User Interface (UI).

NOTE: The UI can be used with gloves.

Cleaning the trophon2 Touch Screen

The screen can be wiped clean using a soft, nonabrasive and lint-free cloth or wipe.

B2.6 Basic Settings

The available settings are:

Sleep timer: The default inactive period until the trophon2 enters sleep mode is two hours, which can be adjusted in this setting.

Language

Date

Time: For record accuracy, user must set the timezone in this setting in addition to the time.

Printer label

Daily timers: Configure a preset daily wake and sleep time

Brightness and Alarm

Network

To access:

1. Select **Menu → Settings**
2. Select the required setting and follow the on screen prompts.

B2.7 AcuTrace™

AcuTrace is an automated HLD traceability system that incorporates Radio Frequency Identification (RFID). AcuTrace enabled accessories and consumables contain RFID chips that store digital information and can be read by the in-built AcuTrace reader.

The trophon2 AcuTrace Reader is located on the device as shown:

AcuTrace compatible products that are embedded with an RFID chip can be identified by this symbol:

Scan the symbol on the AcuTrace reader when prompted.

Products that are AcuTrace compatible have this symbol on the outer packaging:

This symbol cannot be read by the AcuTrace reader.

The following products are compatible with AcuTrace:

trophon AcuTrace Operator Card

trophon AcuTrace Operator Card links the HLD cycle to the operator responsible for various aspects of the workflow. The trophon2 screen will prompt the operator to scan their Operator Card against the AcuTrace reader when preparing a probe for disinfection.

trophon AcuTrace Medical Instrument Tag

trophon AcuTrace Medical Instrument Tag stores probe identification information. Attach the trophon AcuTrace Medical Instrument Tag at a distance from the probe that will not interfere with the patient exam or closing of the trophon door. The tag is scanned on the AcuTrace reader when prompted by an on screen message, prior to HLD, to link the HLD cycle to the disinfected probe.

trophon Chemical Indicator and Disinfectant

The AcuTrace compatible trophon Chemical Indicator box and disinfectant cartridges can be scanned by the AcuTrace reader to track batch and expiry dates of these products.

Refer to each product Instruction for Use (IFU) for more information on their AcuTrace capability and how to programme

B2.8 AcuTrace™ Settings

AcuTrace features can be enabled or disabled via the trophon2 menu. Press:

1. **Menu → AcuTrace → Workflow.**
2. Select the required setting for update and follow the on screen prompts.

SECTION B3: trophon AcuTrace PLUS

trophon **AcuTrace PLUS** allows trophon2 to connect to Hospital Information Systems for automated sharing and centralised storage and processing of disinfection data and enables Parametric Release.

B3.1 Activation

trophon **AcuTrace PLUS** will be prompted during the Initial Setup and can be activated through the product trophon **AcuTrace PLUS** Activation Card. Follow the on screen instructions during the initial setup.

To activate after initial setup, select:

Menu → AcuTrace → AcuTrace PLUS

Follow the on screen instructions.

B3.2 Network Parameters Setup

Consult your facilities IT administrator or network parameters expert to complete this setting.

The device supports the TCP/IP protocol. Set the IP address using the Dynamic Host Configuration Protocol (DHCP) or a static configuration, select:

Menu → Settings → Network

When configured using DHCP, the device will be assigned a valid IP address from the router.

To set the IP address with the static configuration, please check prior with a network administrator to ensure that the IP address is not in use and advise the appropriate values for the subnet mask and the default gateway.

PART C – OPERATION

SECTION C1: Loading the Disinfectant Cartridge

The cartridge door opens automatically, DO NOT force open.

A trophon2 disinfectant cartridge needs to be inserted into the trophon2 before a HLD cycle can commence.

Refer to the disinfectant IFU for detailed instructions on how to scan and log, insert, or remove disinfectant cartridges.

Inserting a Disinfection Cartridge

The trophon2 screen will automatically prompt to scan (if AcuTrace is enabled) and insert a new disinfectant cartridge, when required. Follow the on screen instructions and refer to the disinfectant IFU.

SECTION C2: Logging the trophon Chemical Indicators

No other chemical indicators are approved for use in the trophon2.

If AcuTrace is enabled, log the new batch of Chemical Indicators by selecting:

1. **Menu → AcuTrace → Log Chemical Indicators.**
2. Follow the on screen instructions.

SECTION C3: Routine High Level Disinfection Cycle

C3.1 Preparing the Probe

Wear gloves throughout the complete HLD process.

Clean and dry the probe and check clearly for probe defects BEFORE commencing the high level disinfection process, as per the probe manufacturer's instructions. Ensure that the probe has been thoroughly dried after cleaning and no visible debris is present.

C3.2 Inserting the Chemical Indicator

A whole trophon Chemical Indicator must be used for each disinfection cycle and may only be used once. Refer to the Chemical Indicator IFU.

C3.3 Positioning the Probe

1. If AcuTrace is enabled, when prompted, scan the trophon AcuTrace Medical Instrument Tag and Operator Card.
2. Open the chamber door and load probe and Chemical Indicator.
3. Two clamps (see Figure 9) hold the probe cable securely in the chamber.
4. Wear gloves and insert the probe into the trophon2 by gently pulling the probe cable against the cable clamp (see Figure 9a). Then carefully pull the cable upwards until the probe is suspended in the correct location and the probe cable is held by the cable seal (see Figure 9b). Refer to figures 10-12 and note a correctly positioned probe in the trophon2.

Do not pull the probe down when loaded in the cable clamp.



Figure 9

5. Ensure that the probe is correctly positioned in the chamber. The probe must not make contact with the chamber wall and must be positioned at or above the embossed line at the bottom of the chamber.
6. If the probe is disconnected from its connector end, the external portion of the cable and connector can be secured by placing the connector carefully inside the cable tray.

NOTE: Incorrect positioning of the probe may result in:

- Unsuccessful HLD cycles.
- Residual disinfectant remaining on the probe's surface. This may lead to temporary bleaching and / or irritation of the skin if gloves are not worn.
- Damage to the probe may occur if in contact with the chamber wall.

NOTE: Curved probes must be correctly inserted in the trophon2. See Figure 12.

Curved Probes

The Integrated Probe Positioner (IPP) improves the positioning of approved curved probes in the trophon2 chamber (refer to the Validated Probe List).

Position the probe curve towards the **left** of the trophon2 chamber. To engage the IPP, once the probe is inserted into the probe clamp, push the left side of the IPP and slide to the right (see Figure 13a). This should displace the probe cable till the probe is no longer in contact with the chamber wall (Figure 12).

To release the IPP, depress the trigger (see Figure 14a) and slide to the left.


Figure 13


Figure 14

C3.4 Closing the Chamber Door

- The chamber door utilises a two stage closure mechanism. Carefully close the chamber door to the first click and do NOT force shut. The door will automatically close further to seal and lock at the start of the HLD cycle.
- If the door is not properly closed, a screen message will prompt to close the chamber door.

DO NOT force or manually unlock the chamber door during a HLD cycle.

C3.5 Disinfecting the Probe

Once the door is closed, the trophon2 will request confirmation that the probe is clean and dry prior to HLD. The door will lock upon confirmation.

1. Scan the trophon AcuTrace Operator Card when prompted (if AcuTrace enabled) and follow the on screen instructions.

2. If the probe has NOT been cleaned and dried, press *No*. Remove the probe, clean and dry the probe and repeat the above steps from C3.3.
3. Once the probe is cleaned and dried and the probe has been correctly inserted into the trophon2, press the cycle start symbol on the UI to initiate the cycle or *Cancel* to unlock the machine, remove the probe and follow the on screen instructions.
4. The progress of the disinfection cycle is indicated on the screen.
5. The HLD cycle will take 7 minutes to complete.

In the unlikely event that mist is escaping from the chamber, avoid direct contact with the mist and remain at distance from the trophon2 until completion of the operating cycle and the mist is no longer visible. Contact your customer service representative (see Part D – Troubleshooting).

C3.6 Removing the Probe

After completion of a successful HLD cycle, the ultrasound probe and chamber may have surface temperatures up to 45°C/ 113°F and 60°C/ 140°F respectively. Care should be taken not to touch the chamber. The probe will feel warm to touch but is safe for handling and use with gloves on.

1. Wear gloves and follow the on screen instructions.
2. Remove and wipe the probe with a dry, clean, single-use, low lint or lint free cloth. Visually inspect the probe and remove any present disinfectant residue.

NOTE: Take care when removing the probe to ensure minimal contact - avoid touching the probe against the chamber's hot surface and do not compromise the probe's high level disinfected surface before use.

3. Remove the used Chemical Indicator and verify the colour change against the chart on the Chemical Indicator packaging. Discard the used Chemical Indicator.
4. Close the chamber door and follow the on screen instructions to record the Chemical Indicator result.
5. trophon2 will also indicate a successful or failed cycle. If the Chemical Indicator and / or trophon2 indicated a failed cycle, repeat steps from section C3.1.
6. Discard gloves.
7. The probe is now ready for use or storage in a Nanosonics' Clean Ultrasound Probe Cover.

C3.7 Sleep Mode

If the trophon2 is inactive for an extended time period, it will automatically enter the power saving sleep mode, indicated by the on screen symbol:

Touch the symbol to reactivate the trophon2 for use.

To adjust the trigger time for Sleep Mode, select: **Menu → Settings → Sleep timer**

PART D – RECORDS

SECTION D1: Record Options

All disinfection cycles are automatically logged on trophon2.

If AcuTrace is set up, records of operators, probes, Chemical Indicators and disinfectant are also logged. The most recent records can be viewed on screen or the complete records downloaded to USB.

If a trophon Printer is connected to trophon2, a cycle record will be automatically printed when the HLD cycle is completed. Up to 4 records can be printed, however, 2 records are printed by default.

Refer to the trophon Printer User Manual for more information.

To access the records, select:

Menu → Records then choose to view one of the record list options detailed below and follow the on screen instructions.

Last Cycles

View most recent HLD cycles on screen and optionally print the details of a selected cycle.

Disinfectant *

Disinfectant cartridges scanned and used by the trophon2.

Chemical Indicator *

Chemical Indicator batches logged and used by the trophon2.

Probe *

Probes scanned for disinfection in the trophon2.

Operator *

trophon2 Operators logged by scanning their Operator Card.

NOTE: The Operator ID listed (if AcuTrace enabled) is the user responsible for logging the pass or fail Chemical Indicator result via the Operator Card at the end of HLD cycle.

USB Export

Disinfection cycle records can be downloaded onto USB. Insert a USB drive into any of the three ports on the left hand side of the trophon2, then select **Download records to USB** and follow the on screen instructions.

*AcuTrace must be enabled to view this record

PART E –MAINTENANCE AND ROUTINE CARE

SECTION E1: Preventative Maintenance Service

Your trophon2 requires annual preventative maintenance and a major service at every 5,000 cycles.

Once the service due reminder screen message is displayed, please contact your customer service representative to arrange a service. Contact details, if not provided by your distributor can be found on the Nanosonics website.

The service options found: **Menu → Maintenance** are used to troubleshoot device malfunctioning and should only be used when advised by service personnel.

SECTION E2: Purge Cycle

The purge cycle removes and converts any remaining disinfectant from the trophon2, into oxygen and water.

E2.1 When to Run a Purge Cycle

- Prior to transporting the trophon2, if it has previously been in use.
- When an on screen message states that your trophon2 requires purging. This will occur upon disinfectant expiry at 30 days after insertion. Follow the on screen instructions to purge.

NOTE: Purging can be deferred until convenient on the day of expiry.

Initiate a purge cycle:

- When the trophon2 detects an error that requires a service call.
- Before lifting or moving the trophon2.
- For troubleshooting purposes when directed by service support only.

E2.2 How to Initiate a Purge Cycle

NOTE: Once the purge cycle has started it may be paused but it cannot be cancelled. Do not switch the trophon2 off during purging as this will restart the purge cycle. Do NOT attempt to open the chamber or cartridge door during the purge cycle.

To initiate a purge cycle:

1. Wear gloves and ensure that the waste drawer is empty and fully inserted into the trophon2 (see Figure 3 for location).
2. Select: **Menu → Maintenance → Purge**. Follow the on screen instructions.
The purge cycle will typically take less than 30 minutes.
3. When purging is complete, put on gloves, empty the waste drawer and follow the onscreen messages.
4. Add a new disinfectant cartridge (unless transporting), follow the on screen instructions and refer to the disinfectant IFU.

SECTION E3: Regular Cleaning

1. Do NOT submerge the trophon2, or pour liquids over it.
2. Keep the trophon2 level and upright at all times.
3. Keep the power socket completely dry (see Figure 4).
4. For cleaning, when cool, wipe the chamber and the outside surfaces of the trophon2 with a cloth or wipe moistened with a mild, soapy solution, until all surfaces are visibly clean.
5. For disinfecting, wipe all accessible surfaces of the trophon2 with an Isopropanol or Quaternary Ammonium (Quat) Wipe.

SECTION E4: Transporting the trophon2

NOTE: The below is not necessary for transportation of the device using the trophon Cart.

To transport the trophon2:

- You must purge the disinfectant (do not insert a new disinfectant cartridge until trophon2 is relocated), switch off the trophon2 power switch and unplug from the mains.
- Keep the trophon2 upright at ALL times.
- Pack using only Nanosonics approved packaging.

SECTION E5: Disposal of trophon2

Contact your Nanosonics representative to advise the appropriate collection point for the recycling of electrical and electronic equipment.

PART F –TROUBLESHOOTING

SECTION F1: Incomplete or Failed Cycles

This section describes the most common causes for a cycle malfunction and the advised steps to action.

F1.1 Mains Power Failure

If the mains power supply to the trophon2 is lost whilst in operation, the current cycle will not complete.

- Once the power is restored, follow the on screen messages to remove the probe safely from the trophon2.
- Discard the used Chemical Indicator and replace with a new one.
- Repeat the HLD cycle.
- If power cannot be restored and the probe is urgently required, follow Section F2.

F1.2 Cycle Fault

If an error occurs during or at the end of a cycle, a cycle fault will be detected. Follow the on screen messages to fix and repeat the HLD cycle.

In case of a repeated fault or serious malfunction, note the screen error message and contact your customer service representative. Do NOT attempt to use the trophon2 or the probe.

SECTION F2: Manual Door Lock Override

Use **ONLY** when the probe is locked in the chamber and must be **urgently** retrieved.

Chamber surfaces may be hot and disinfectant present.

Gloves must be worn to avoid contact with disinfectant.

Do NOT manually open the door during a HLD cycle. Disinfectant mist will be present and contact must be avoided.

1. Ensure the trophon2 is powered off.
2. Open the waste drawer, remove any liquid and obtain the key (Figure 15).
3. Locate the Manual Door Lock Opening Mechanism Cover on the right side of the trophon2 (Figure 16). Align the 4 key notches with the indents on the cover and turn the key **ANTI CLOCKWISE** to unscrew.
4. Once the Manual Door Lock Opening Mechanism Cover is removed, align the 4 notches with the grooves of the Manual Door Lock Opening Mechanism inside. Push and turn **CLOCKWISE** to 90 degrees to unlock the chamber door (Figure 17).

The probe is NOT DISINFECTED and CANNOT be used until it has completed a successful HLD cycle or high level disinfected by an alternative method.



Figure 15

Figure 16

SECTION F3: Diagnostics

Use this table to diagnose and resolve basic problems. If a probe requires urgent retrieval from trophon2, follow Section F2.

If the problem persists, contact your trophon2 representative.

Symptom	Check / action the following:
1. There is no power to the trophon2. 2. The screen is blank.	<ul style="list-style-type: none"> The trophon2 is completely plugged in and switched ON at both the trophon2 power switch and at the wall outlet. The correct power cable for your region is used.
3. The chamber door will not open.	<ul style="list-style-type: none"> There is power to the trophon2. There is no HLD, warm up or purge cycle running.
4. The chamber door will not close or lock.	<ul style="list-style-type: none"> The probe is loaded correctly. There is nothing obstructing the chamber door or locking mechanism
5. The chamber door is open and locked.	<ul style="list-style-type: none"> Power the trophon2 OFF and back ON using the power switch (see Figure 4).The trophon2 should then automatically unlock the chamber door. If the above are not successful, power OFF the trophon2 and follow section F2.
6. The cartridge door does not open.	<ul style="list-style-type: none"> There is power to the trophon2. A HLD cycle is not in progress. The waste drawer is empty and fully inserted. There is nothing obstructing the cartridge door. <p>Note: the cartridge door is automatic, and should not be forced open.</p>
7. The cartridge door will not close.	<ul style="list-style-type: none"> Correct cartridge type has been inserted. Cartridge is correctly positioned. Cartridge lid has been removed.
8. The probe will not sit correctly in the chamber.	<ul style="list-style-type: none"> The probe is validated for use in the trophon2 – see section A1.2. The probe is loaded correctly. The IPP is correctly engaged – see section C3.3.
9. The cycle will not start.	<ul style="list-style-type: none"> The probe is loaded correctly. The chamber door is closed. Confirm probe is clean & dry before pressing start. All on screen instructions have been followed correctly.

10. Liquid is leaking from the trophon2.	 WARNING: Any fluid leaking from the trophon2 may contain hydrogen peroxide. If liquid or mist is seen coming from the trophon2 at any time: <ol style="list-style-type: none"> 1. Do not come into contact with the mist or liquid. 2. Wear appropriate PPE. 3. Ensure area is well ventilated. 4. Allow the trophon2 to complete the cycle. 5. Turn off the trophon2 and remove the power cord. 6. Contact your customer service representative. 7. Consult the Safety Data Sheet.
11. The trophon2 is failing multiple cycles.	<ul style="list-style-type: none"> • Record any error codes and colour of the Chemical Indicator and contact your customer service representative.
12. The trophon2 takes too long to warm up between cycles	<ul style="list-style-type: none"> • Ensure that the probe is removed promptly at the end of each cycle.

PART G – SERVICE REQUIREMENT AND WARRANTY PROVISION

Contact your trophon2 representative if you have any questions about:

- The trophon2, consumables or accessories
- The warranty

Each trophon2 has a comprehensive warranty against defects in material and workmanship for 12 months from the Date of Purchase. The specific warranty terms and conditions are defined in appendix 2 of this manual. Please be aware of the exclusions.

To ensure the safety and efficacy of your HLD operations, trophon2 requires a preventative service every 12 months and a major service at every 5000 cycles.

Service Schedule: When trophon2 is due for service, a service due message will be displayed on the screen to prompt service arrangements. The message will be indicated prior to commencing a HLD cycle at weekly intervals until a service is performed. Service due information can also be accessed by selecting:

Menu → Maintenance → Service data.

Nanosonics has made the service provision available to customers through either our direct service or our service partners including local distributors who have been trained and authorised to service trophon2. Only authorised customer service or suitably trained personnel should service the trophon2 with genuine parts supplied from Nanosonics.

Modifying the trophon2 without authorisation will void your warranty.

APPENDIX 1: TROPHON2 TECHNICAL SPECIFICATIONS N05000-2

N05000-2 Electrical Specification	Rated mains input voltage: 230V AC Rated mains input current: 6Amp, 50/60Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using the power cable supplied with the trophon2.
Data port	Ethernet connector RJ45 USB Port: Type A
Environmental Specification	Operating temperature range: 17 to 27°C / 62.6°F to 80.6°F
Storage and Transport Conditions	Temperature range: -20°C to +60°C / -4°F to +140°F
Physical Characteristics	Weight of trophon2: Unpacked: 22kg (48.5lbs) Packed: 25kg (55lbs) Dimensions of trophon2: 535mm high x 360mm wide x 317mm depth (21in high x 14.2in wide x 12.5in depth)
Electromagnetic Compliance	The trophon2 has been tested and found to comply with the limits for emission (electromagnetic Interference) pursuant to EN61326-1:2013 (CISPR 11 Group 1 Class B limits)

APPENDIX 2: PRODUCT WARRANTY TERMS AND CONDITIONS

Terms

This warranty is given by Nanosonics Limited ABN 11 095 076 896 of 14 Mars Road, Lane Cove, NSW 2066 Australia (**Nanosonics**).

Nanosonics warrants to the customer that the trophon2 is free from defects in material and workmanship that materially affect its function under normal use and service for a period of 12 months commencing upon the Date of Purchase (**warranty period**).

Exclusions

This warranty does not apply in the following circumstances (regardless of how those circumstances arise):

- a. where the trophon2 has not been used, handled, installed, stored, cleaned and serviced in accordance with the relevant user manual or other written instructions issued by Nanosonics (including where used in temperature or other external conditions exceeding those set out in the product specification, or serviced by persons other than Nanosonics' approved service personnel);
- b. where modifications have been made to the trophon2, other than by Nanosonics or its authorised service providers;
- c. where unauthorised consumables, accessories or other chemicals or items have been used with the trophon2;
- d. where authorised consumables, accessories or other chemicals or items have been used inappropriately or incorrectly with the trophon2;
- e. where the trophon2 is used in conjunction with other equipment or products (other than multiple use ultrasound probes as described in the user manual), without Nanosonics' prior written consent;
- f. where the trophon2 has been damaged due to external or environmental causes of any kind (including factors such as voltage fluctuations, excess voltage or power failure);
- g. where the trophon2 has been damaged as a direct or indirect result of any malicious or negligent act or omission by any person (other than Nanosonics or its authorised service providers);
- h. where the defect does not materially affect the function of the trophon2 (for example scratches or marks on the external surface of the trophon2); or
- i. where the serial number or product label has been removed, changed, deleted or made unrecognizable, or if the number or label is no longer clearly distinguishable for other reasons beyond Nanosonics' control and therefore it is not possible to conclusively identify the product; or
- j. where you have not followed a reasonable instruction of Nanosonics however communicated to you

This warranty applies to the trophon2 and defective parts only; the warranty does not cover the replacement of used disinfectant cartridges or of parts which need periodic replacement during the life of the product as a result of the ordinary use made of them.

How to make a claim

Please contact your trophon2 customer service representative with any queries regarding this warranty or post warranty repairs. If you wish to make a warranty claim, please contact your trophon2 customer service representative.

Nanosonics will make arrangements for and bear the cost of the collection of your trophon2. You will be responsible for uninstalling, reinstalling and recommissioning the trophon2, regardless of whether or not it is found to be defective. If Nanosonics finds on examination that the trophon2 is defective in materials and workmanship and is within the warranty period, then we will repair or replace the defective trophon2 at our discretion. Nanosonics will bear the reasonable cost of return delivery of the repaired trophon2, or replacement trophon2, to you at a time stipulated by Nanosonics. If we are unable to repair or replace the trophon2 for any reason, we will discuss with you an appropriate solution including upgrading you to a newer model or refunding the purchase price.

If Nanosonics finds on examination that the trophon2 is **not** defective in materials and workmanship in Nanosonics sole opinion, or if you are not entitled to the benefit of this warranty (for example, if any of the above exclusions apply, or the claim was not made within the warranty period), then Nanosonics may require you to bear the cost of returning the trophon2 to you, and the costs of any repairs to the trophon2, or replacement trophon2, requested by you.

You are responsible for backing up all data on the trophon2 if it is being repaired, and acknowledge that repair of the trophon2 may result in the loss of user-generated data stored on the trophon2. You accept that Nanosonics is not responsible for any loss of data.

Goods presented for repair may be replaced by refurbished goods of the same type rather than being repaired. Refurbished parts may be used to repair the goods.

Australia: *The following statement only applies if you are a ‘consumer’ for the purpose of the Australian Consumer Law at Schedule 2 of the Competition and Consumer Act 2010.*

Our goods come with guarantees that cannot be excluded under the *Australian Consumer Law*. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. The benefits under this product warranty are in addition to other rights and remedies you may have under law in relation to our goods.

New Zealand: *The following statement only applies if you are a ‘consumer’ for the purpose of New Zealand’s Consumer Guarantees Act 1993.*

Our goods come with guarantees that cannot be excluded under the *Consumer Guarantees Act 1993*. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

United States: *The following statement only applies to purchasers of the trophon2 in the United States.*

This warranty is a limited warranty, is the sole and exclusive warranty applicable to the product(s) described herein and is made in lieu of all other warranties, express or implied, including without limitation warranties of merchantability or fitness for a particular purpose.

Other Countries: You may have statutory rights in relation to the trophon2 and these are not affected by this warranty.

APPENDIX 3: MicroDoc LICENCE AGREEMENT

By purchasing and/or installing and/or using the trophon2, you accept and agree to be bound by the terms of the below End User License Agreement (**EULA**) relating to the use of MicroDoc Oracle Java SE Embedded 8 software (the **Software**) utilised in the trophon2. If you do not agree to all of the terms of the EULA, you must not install and/or use the trophon2.

The capitalised terms below have the same meaning as appear in the EULA entered into between Nanosonics and MicroDoc GmbH for the Software (a copy of which is available on request).

- The user is permitted to the use of the MicroDoc Deliverables only according to the *Product Registration Form* (which may be provided on request);
- The user is prohibited from redistributing the *MicroDoc Deliverable*;
- The user is prohibited from unauthorized duplication of the *MicroDoc Deliverables* except for backup or archival purposes;
- The user is prohibited from causing or permitting the translation, de-compiling, disassembly, reverse engineering, unbundling or extraction of the *MicroDoc Deliverables*;
- The user is prohibited from removal or altering of any proprietary notices, labels or marks in or on the *MicroDoc Deliverables*;

- The user is prohibited from export of the *MicroDoc Deliverables* in contravention of applicable export control laws;
- The *End User* acquires no right, title or interest in or to the *MicroDoc Deliverables*;
- The *End User* may only copy the *MicroDoc Deliverables* as necessary to use the *MicroDoc Deliverables* according to the *Product Registration Form*, to follow normal archiving practices, and shall use reasonable efforts to ensure that all copies of the *MicroDoc Deliverables* bear any notice contained on the original *MicroDoc Deliverables*;
- *LICENSEE'S* licensors are intended third party beneficiaries of all terms and conditions that apply to the *MicroDoc Deliverables*, including all warranty and liability limitations and any right of indemnification, and a list of such beneficiaries will be provided upon written request.
- The *End User* is prohibited from creating, modifying, or changing the behavior of classes, interfaces, or subpackages that are in any way identified as "Java", "Javax", "Sun" or similar convention as specified by Oracle in any naming convention designation;
- The *End User* acknowledges that Oracle owns the Java trademark and all Java-related trademarks, logos and icons including the Coffee Cup and Duke ("Java Marks") and agree to: (a) comply with the Java Trademark Guidelines at <http://www.oracle.com/us/legal/third-partytrademarks/index.html>; (b) not do anything harmful to or inconsistent with Oracle's rights in the Java Marks; and (c) assist Oracle in protecting those rights, including assigning to Oracle any rights acquired by the *End User* in any Java Mark;
- The *LICENSEE* hereby notifies the *End User* that the Programs may contain source code that, unless expressly licensed for other purposes, is provided solely for reference purposes pursuant to the terms of the *End User* license agreement;
- The *LICENSEE* hereby notifies the *End User* that additional copyright notices and license terms applicable to portions of the Programs are set forth in the THIRDPARTYLICENSEREADME.txt file;
- The *LICENSEE* hereby notifies the *End User* that use of the Commercial Features for any commercial or production purpose requires a separate license from Oracle. "Commercial Features" means those features identified in Table 1-1 (Commercial Features In Java SE Product Editions) of the Program Documentation accessible at <http://www.oracle.com/technetwork/indexes/documentation/index.html>
- The *END USER* may only distribute the *MicroDoc Deliverables* only as incorporated in and as an integral part of their product according to the *Product Registration Form*.

APPENDIX 4: TROPHON2 ACCESSORIES*

The additional accessories below are designed for use with the trophon2 and available for purchase:

trophon AcuTrace Operator Card

Operator Card can be scanned on the trophon2 to link the HLD cycle to the operator responsible for various aspects of the workflow.

trophon AcuTrace Medical Instrument Tag

Medical Instrument Tags are attached to the probes and scanned at the start of a HLD cycle to link HLD to the disinfected probe.

trophon AcuTrace Plus Activation Card

An add on package that allows trophon2 to connect to Hospital Information Systems for automated sharing and centralised storage and processing of disinfection data. This package also enables the Parametric Release functionality of the trophon2.

trophon Wall Mount 2

The wall mount allows appropriate wall attachment of trophon2.

trophon Cart

The cart permits mobility of the trophon2.

trophon Clean Ultrasound Probe Cover

Disposable plastic covers that protect and maintain disinfected probes between uses.

trophon Printer & Label Roll

The trophon2 compatible printer uses the label roll to print disinfection records.

trophon Printer Wall Mount

Suitably attaches the trophon Printer to the wall.

trophon Printer Cart Mount

Allows attachment of the trophon Printer to the trophon Cart.

trophon Companion Drying Wipes

Multi-purpose drying wipe

Only use trophon approved accessories or the trophon2 may be ineffective.

GLOSSARY**AcuTrace™**

The RFID technology used by trophon2 and associated products as part of the disinfection traceability system.

Cable clamp

Mechanism at the top of the trophon chamber to grip and hold the probe cable during HLD.

Cartridge

The disinfectant product container that is inserted into the trophon2.

Chemical Indicator

A consumable that detects and indicates the correct concentration of disinfectant during a cycle to ensure HLD.

Disinfectant

The consumable liquid in cartridge used by trophon2 responsible for HLD.

Cable Seal

A sealing mechanism at the top of the trophon2 chamber (below the cable clamp) to prevent disinfectant leaking from the chamber.

High Level Disinfection (HLD)

The treatment of semi- critical medical devices with disinfectant to inactivate all microorganisms, excluding spores.

Instruction for Use (IFU)

Instruction guide for recommended use of the product.

Integrated Probe Positioner (IPP)

Mechanism inside the trophon2 chamber to assist with the correct positioning of probes during HLD.

Minimum Effective Concentration (MEC)

The minimum concentration of the disinfectant required for HLD.

Purge Cycle

The process of removing all disinfectant from the trophon2 system.

RFID

Radio Frequency Identification

Semi-critical Devices

Medical devices that contact mucous membranes and non-intact skin but do not penetrate sterile tissue.

Safety Data Sheet

An outline of the potential health risks posed by a hazardous chemical and the safe working procedures required to mitigate those risks.

User Interface (UI)

The human controlled interaction with software or machine.

Validated Probe List

A list of probes that have been tested and approved for use in the trophon2 by Nanosonics, in conjunction with the probe Original Equipment Manufacturers (OMEs).

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Nanosonics, Inc.
Nanosonics Europe GmbH