trophon[®]

Chemical Indicator

 Nanosonics Limited
14 Mars Road, Lane Cove
NSW 2063, Australia
+61 2 8063 1600
info@nanosonics.com.au
www.nanosonics.com.au

USA Contact:
Nanosonics, Inc.
7205 E. 87th Street
Indianapolis, Indiana 46256
1-844-876-7466
support@nanosonics.us
www.nanosonics.us

 Nanosonics Europe GmbH
Poppenbuetteiler Bogen 66
22399 Hamburg - Germany
+49 40 46856885
www.nanosonics.eu

 Nanosonics Limited UK
Unit 2 Linfit Court, Colliers Way,
Clayton West, Huddersfield,
HD8 9WL
+44 1484 860581
ukinfo@nanosonics.co.uk
www.nanosonics.co.uk

E Cat: E8350MB

REF N05003

EN - Symbol Glossary / DE - Glossar der Symbole / FR - Glossaire des symboles / NL - Verklarende woordenlijst van symbolen / SV - Symbolförklaring / JP- 記号の説明 / ES - Glosario de símbolos / IT - Glossario dei simboli / PT - Glossário de símbolos / KR - 기호 설명 / RU - Перечень символов / NO - Symbolordliste / DK - Symbolordliste / FI - Symbolien selitykset

EN	DE	FR	NL	SV	JP	ES	IT	PT	KR	RU	NO	DK	FI
EN - Date of Manufacture DE - Herstellungsdatum FR - Date de fabrication NL - Fabricagedatum SV - Tillverkningsdatum JP - 製造日 ES - Fecha de fabricación IT - Data di produzione PT - Data de fabrico KR - 제조일 RU - Дата изготовления NO - Produktionsdato DK - Fabricationsdato FI - Valmistusajankohhta	EN - Manufacturer DE - Hersteller FR - Fabricant NL - Fabrikant SV - Tillverkare JP - 製造者 ES - Fabricante IT - Fabricante PT - Fabricante KR - 제조업체 RU - Производител NO - Produsent DK - Fabrikant FI - Valmistaja	EN - Product Number DE - Artikelnummer FR - Code produit NL - Productnummer SV - Produktnummer JP - 製品番号 ES - Número de producto IT - Numero prodotto PT - Referência KR - 제품 번호 RU - № изделия NO - Produktnummer DK - Produktnummer FI - Tuotenumero	EN - Batch No. DE - Chargennr. FR - N° de lot NL - Partijnr. SV - Batchnr. JP - バッチ番号 ES - N.º de lote IT - N. partita PT - N.º de lote KR - 배치 번호 RU - № партии NO - Batchnummer DK - Batch nr. FI - Eränumero										
 EN - Waste Electrical and Electronic Equipment DE - Abfall elektrische Ausrüstung FR - Déchets d'équipements électriques et électroniques NL - Afgedankte elektrische en elektronische apparatuur SV - Elektriskt och elektroniskt avfall JP - 電気・電子機器廃棄物 ES - Eliminación de equipos eléctricos y electrónicos IT - Rifiuti di apparecchiature elettriche o elettroniche PT - Resíduos de Equipamentos Elétricos e Eletrónicos KR - 전기 및 전자 장비 폐기물 RU - Отходы электрического и электронного оборудования NO - Avfall av elektrisk og elektronisk utstyr DK - Affald af elektrisk og elektronisk udstyr FI - Sähkö- tai elektroniikkajätteen häilytys	 EN - Keep Out of Direct Sunlight DE - Vor direktem Sonnenlicht schützen FR - Conserver à l'abri de la lumière directe du soleil NL - Uit de buurt van direct zonlicht houden SV - skydd från solen JP - 直射日光を避けること ES - Mantener apartado de la luz directa del sol IT - Tenere al riparo dalla luce solare diretta PT - Manter afastado da luz solar KR - 직사광선 주의 RU - Не подвергать воздействию прямых солнечных лучей NO - Oppbevares borte fra direkte sollys DK - Hold væk fra sollys FI - Säilytä pois auringonvalosta												

EN	DE	FR	NL	SV	JP	ES	IT	PT	KR	RU	NO	DK	FI
EN - European Authorised Representative DE - Bevollmächtigter FR - Mandataire européen NL - Geautoriseerde vertegenwoordiger in Europa SV - Europeisk auktoriserad representant JP - 欧州代表 ES - Representante europeo autorizado IT - Mandatario per l'Europa PT - Mandatário europeu KR - 유럽대리인 RU - Уполномоченный представитель ЕС NO - Europeisk autorisert representant DK - Europæisk autoriseret repræsentant FI - Valtuutettu edustaja Euroopassa	EN - AcuTrace® RFID Zone FR - Zone RFID AcuTrace® NL - Zone RFID AcuTrace® SV - Zonen RFID AcuTrace® JP - RFIDゾーン ES - Zona RFID AcuTrace® IT - Zona RFID AcuTrace® PT - Zona RFID de AcuTrace® KR - RFID 구역 RU - Зона RFID AcuTrace® NO - AcuTrace® RFID-zone DK - AcuTrace® RFID-zone FI - AcuTrace® RFID -alue	EN - EU Importer FR - Importateur UE FR - Zone d'expiration NL - Importeur in EU SV - Importör i EU JP - EU輸入者 ES - Importador de la UE IT - Importatore UE PT - Importador da UE KR - EU 수입업체 RU - Импортёр ЕС NO - EU-importør DK - EU-importør FI - Maahantuaja Euroopan Unionin alueella	EN - Expires DE - Verfalldatum FR - Date d'expiration NL - Utløper SV - Utgåår JP - 使用期限 ES - Caducida IT - Scadenza PT - Prazo de validade KR - 유효 기한 RU - Истекает NO - Går ut på dato DK - Vældstid FI - Vanhee										

			
EN - Single use DE - Für den einmaligen Gebrauch FR - À usage unique NL - Voor eenmalig gebruik SV - För engångsbruk JP - 単回使用 ES - De un solo uso IT - Monouso PT - Não reutilizável KR - 일회용 RU - Для одноразового использования NO - Engangsbruk DK - Enkelt brug FI - Kertaikäyttöinen	EN - Instructions for Use DE - Gebrauchsanweisung FR - Qualification de l'Installation NL - Gebruiksaanwijzing SV - Bruksanvisning JP - 使用法 ES - Calificación de la instalación IT - Qualifica per l'installazione PT - Qualificação para instalação KR - 설치 자격 RU - Установочная квалификация NO - Installasjonskvalifisering DK - Kvalificatie van installatie FI - Asennus, ammattitaitovaatimus	EN - Storage temperature range between 15°C– 30°C(59°F – 86°F) DE - emperaturbereich für die Aufbewahrung 15–30 °C FR - température de conservation : entre 15 °C et 30 °C NL - Opslagtemperatuur tussen 15°C–30°C SV - emperaturmråde för förvaring 15–30 °C JP - 保管温度範囲 15℃～30℃ ES - Intervalo de temperaturas dealmacenamiento: 15°C– 30°C (59°F – 86°F) IT - Intervallo da temperatura tra 15 e 30°C PT - Intervalo da temperatura de armazenamento entre 15 °C e 30 °C KR - 보관 온도 범위: 15℃ – 30℃ RU - температура хранения 15°C–30°C NO - Lagringstemperaturområdet 15 °C – 30 °C DK - Opbevaringstemperatur 15 °C – 30 °C FI - Säilytyslämpötila-alue 15–30 °C	

EN Instructions for Use
Der trophon[®] Chemical Indicator is a qualitative chemical indicator designed to monitor the minimum effective concentration (MEC) of the disinfectant during the trophon[®]EPR / trophon[®]2 high level disinfection (HLD) process. Only use this product with trophon EPR / trophon2. Please read the trophon EPR / trophon2 User Manual for important information and warnings prior to using the trophon EPR / trophon2 and Chemical Indicator.

Important Information
Use MEC color with enclosed Lot of Chemical Indicators only. Frequency of Use: Use one Chemical Indicator with each trophon EPR / trophon2 HLD cycle. MEC: For the Chemical Indicator to be effective, it must pass the MEC. Refer to the Color Assessment Chart on the packaging for the pass and fail criteria. Pass Criteria: Each HLD cycle is ONLY successful when BOTH: • The color of the Chemical Indicator falls within the pass criteria specified on the Color Assessment Chart (Refer to packaging) AND • The trophon EPR / trophon2 indicates that the cycle is complete. Storage: • Store at room temperature 15°C – 30°C (59°F – 86°F). • Store in a dry and clean environment, out of direct heat and UV light. • Do not store near or expose to chemicals such as sterilising agents, acids, bases, bleaches and other disinfectants.
AcuTrace[®] : Refer to the trophon EPR / trophon2 User Manual and on screen instructions to enable the AcuTrace Technology. Quantity: 300 per pack Shelf Life/Expiry: • Refer to packaging for expiration date. • Discard any unused Chemical Indicators after the expiration date. General Precautions: • Not intended for use with any other HLD process. • Gloves must be worn when removing the Chemical Indicator from the chamber. • Discard used/exposed Chemical Indicator – do not store for later reference. • Do not use Chemical Indicator after expiration date as shown on packaging. • Shipping information: Not Restricted. • The trophon Chemical Indicator is NOT a hazardous product according to the OSHA Hazard Communication Standard 29 CFR 1910.1200. • Discard Chemical Indicator packaging once empty.

- Instructions for Use**
- If AcuTrace is enabled, scan the AcuTrace RFID zone on the Chemical Indicator packaging when a new pack is used (figure 1).
 - Remove a trophon Chemical Indicator by pressing out from the perforated sheet (figure 2) and re-dose the packaging.
 - Ensure the Chemical Indicator is **RED** in color.
 - Place the Chemical Indicator (**red side facing up**) into the holder located at the bottom of the chamber door (figure 3).
 - After loading the ultrasound probe, close the chamber door and follow the onscreen messages and **trophon EPR / trophon2 User Manual** to initiate the HLD cycle.
 - When the HLD cycle is complete, wear gloves and remove the Chemical Indicator from the trophon EPR / trophon2.
 - Visually inspect the Chemical Indicator for color change and verify this against the **Color Assessment Chart** (see packaging). The Chemical Indicator result and trophon EPR / trophon2 screen message should be interpreted in conjunction with each other. Refer to Table 1 for further instructions
 - Discard the used Chemical Indicator.

Following exposure to the trophon EPR / trophon2 HLD process, the trophon Chemical Indicator will change color, indicating a passed or failed cycle:

Table 1. Color Assessment Instructions (Refer to chart on packaging)

Indicator Color	Screen Message	Process Pass/ Fail	Instructions
More yellow than MEC	HLD Cycle Complete	✓ Pass	Pass: - Wear gloves and discard used Chemical Indicator. - Wipe the probe with a dry, clean, single-use, low lint cloth. - The probe is now ready for use or storage.
More yellow than MEC	HLD Cycle Failed	✗ Fail	Fail: - Wear gloves and discard used Chemical Indicator. - Insert new Chemical Indicator. - Repeat the trophon EPR / trophon2 HLD process.
MEC or more red	HLD Cycle Complete	✗ Fail	
MEC or more red	HLD Cycle Failed	✗ Fail	

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