

Customer Information Notice: CIN179

CIN179 trophon®3 compatibility with trophon consumables and accessories

Nanosonics has released the trophon3 device which shares the same core technology of 'sonically activated' hydrogen peroxide of trophon EPR and trophon2 devices.

The trophon consumables and accessories listed in the table below are compatible and useable with trophon3 devices. Nanosonics is updating the product labelling and Instructions for Use of these products to reflect the compatibility with trophon3 devices.

Note that not all products and product codes are available in all countries. Refer to the *Consumables and Accessories* brochure available in your local Customer Centre portal at www.nanosonics.com.

Product Name	Product Code	Quantity
trophon NanoNebulant	N05001	1 Box of 6 Bottles
trophon Sonex-HL	N05002	1 Box of 6 Bottles
trophon Chemical Indicator	N05003	1 Box of 300 Chemical Indicators
trophon AcuTrace Operator Card	N05006	1 Box of 10 Cards
trophon AcuTrace Medical Instrument Tag	N05007	1 Box of 10 Tags
trophon Printer	N00048 N00048-ANZ N00048-NNA	1 Printer
trophon Printer Wall Mount	N00105	1 Mount
trophon Printer Label Roll	N00049-ROW	1 Roll of 600 Labels*

trophon Logbook	N00098	1 Box of 5 Logbooks
trophon Wall Mount 2	N05009	1 Mount
trophon Cart	N05010	1 Box of 1 Cart
trophon Printer Cart Mount 2	N05013	1 Mount
trophon Cart UPS Bracket	N00110	1 Bracket
trophon Wireless Ultrasound Probe Holder (WUPH-A)	N05030 N05030-US	1 Unit
trophon Companion Drying Wipes	N00113	1 Box of 6 Cannisters x 100 Wipes
trophon Companion Wipes Wall Mount	N00114	1 Mount
Clean Ultrasound Probe Covers	N00102	1 Box of 100 Covers

*1 Roll of 350 labels for Parametric Release

Should you have any questions related to what products are compatible with trophon3 devices, please contact your local Nanosonics representative.