

trophon^{®2}

Simply Smarter

High Level Disinfection for Ultrasound Probes



Consumables and Accessories

Now with AcuTrace™ (RFID) technology to enhance digital record keeping and deliver cost effective, easy and reliable workflows for ultrasound probe high level disinfection.



nanosonics
Infection Prevention. For Life.

trophon[®]2

Simply Smarter

High Level Disinfection for Ultrasound Probes



High Level Disinfection Device



trophon2 (N05000-ANZ)

The latest innovation in ultrasound probe high level disinfection (HLD) featuring an enhanced design, simple and fast workflows, and all new AcuTrace™ for digital record keeping and seamless integration with your hospital IT system.



NEW
Shelf life now extended
to 16 months



Consumables



NanoNebulant[®] (N05001)

NanoNebulant is a proprietary disinfectant liquid with 35% hydrogen peroxide chemistry to achieve effective HLD of the entire ultrasound probe including the shaft and handle. Incorporates AcuTrace™ technology to store batch and expiry details directly into trophon2.



trophon Chemical Indicator (N05003)

The trophon Chemical Indicator responds to a combination of specified time, temperature, disinfectant dosage and disinfectant concentration. These critical parameters need to be met for the chemical indicator to change to a pass colour. This qualitative colour change provides independent confirmation of success of each disinfection cycle. Indicator packaging now incorporates AcuTrace™ technology to record and store batch and expiry details directly into trophon2.

Traceability Accessories



trophon AcuTrace™ Operator Card (N05006)

Incorporates AcuTrace™ technology to digitally capture operator interaction details, stored directly on your trophon2.



trophon AcuTrace™ Medical Instrument Tag (N05007)

Incorporates AcuTrace™ technology to digitally capture ultrasound probe details, stored directly on your trophon2.

trophon Printer (N00048)

The trophon Printer is a fast, easy to use traceability solution that helps to link the probe and the HLD cycle with the patient.

trophon Printer USB Cable (N05011)

Cable for connecting the trophon Printer to trophon2. Note: only required if connecting trophon2 to an older trophon printer. All trophon printers now include a USB cable.

trophon Printer Label Roll (N00049-ROW)

High quality printer labels with a minimum 7-year life span, appropriate for medical records.

trophon Logbook (N00098)

Makes recording important disinfection steps easy. Compatible with trophon Printer Labels.



Point of Care Set-up Accessories

trophon Wall Mount 2 (N05009)

A convenient option for installing trophon2 where there are space constraints.

trophon Printer Wall Mount (N00105)

Custom designed for secure, horizontal mounting of the trophon Printer to a wall. Includes a cable tidy area for convenience.

trophon Cart (N05010)

Makes trophon fully mobile for convenient point of care use and adaptability where there are space constraints.

trophon Printer Cart Mount 2 (N05013)

Designed for securely attaching the trophon Printer to the trophon Cart in an easily accessible, convenient location.



Drying Wipes

trophon Companion Drying Wipes (N00113)

Used to dry probes before and after HLD. Unlike many wipes that just move moisture around, trophon Companion Drying Wipes are highly absorbent to leave surfaces completely dry.

trophon Companion Wipes Wall Mount (N00114)

For wall mounting trophon Companion Cleaning and Drying Wipes for easier access and dispensing.



Ultrasound Probe Storage Cover

Clean Ultrasound Probe Cover (N00102)

Custom designed to protect endocavity and surface ultrasound probes from recontamination (through handling or the environment) after HLD and before they are used on the next patient.



IT System Integration Option

trophon AcuTrace™ PLUS Activation Card (N05008)

AcuTrace PLUS™ makes the disinfection records stored on trophon2 accessible to other software applications over a network. Configure your IT systems* to centrally backup all trophon devices on the network, and link individual disinfection cycles with patient records.



*Important note: all connectivity, configuration and integration with customer IT systems is the responsibility of the customer.

trophon2 System Product Descriptions

Product Name	Product Code	Quantity
trophon2	N05000-ANZ	1 Device
NanoNebulant [®]	N05001	1 Box of 6 Bottles**
trophon Chemical Indicator	N05003	1 Box of 300 CIs
trophon AcuTrace [™] Operator Card	N05006	1 Box of 10 Cards
trophon AcuTrace [™] Medical Instrument Tag	N05007	1 Box of 10 Tags
trophon Printer	N00048	1 Printer
trophon Printer USB Cable	N05011	1 Cable
trophon Printer Label Roll	N00049-ROW	1 Roll of 600 Labels [#]
trophon Logbook	N00098	1 Box of 5 Logbooks
trophon Wall Mount 2	N05009	1 Mount
trophon Printer Wall Mount	N00105	1 Mount
trophon Cart	N05010	1 Box of 1 Cart
trophon Printer Cart Mount 2	N05013	1 Mount
trophon Companion Drying Wipes	N00113	1 Box of 6 Can x 100 Wipes
trophon Companion Wipes Wall Mount	N00114	1 Mount
Clean Ultrasound Probe Covers	N00102	1 Box of 100 Covers
trophon AcuTrace [™] PLUS Activation Card	N05008	1 Card

** Minimum order quantity: 2 boxes.

[#] Parametric Release Label: 1 Roll of 350 labels.



The latest innovation in ultrasound probe HLD

featuring an enhanced design, simple and fast workflows, and all new AcuTrace[™] for digital record keeping and seamless integration with your hospital IT system.